

CHANGING THE SETTING FOR DIFFERENT SPECIES

The TONOVET tonometer has measuring settings for dog/cat and horse. The dog and cat use the same setting; displayed as “do”. The horse setting is displayed as “ho”. The third setting “P” is undefined, reserved for future use. The default setting is “do”.

Changing the setting:

1. Press the right or left selector button until the display shows “End”.
2. Press **first** the measurement button, **keep it pressed down and immediately press right selector button**. *Note that buttons have to be pressed down at the same time.*
The display shows “CAL” (or software version).
3. Press the measurement button.
4. Toggle between settings using the selector button.
5. Choose the setting and activate by pressing the measurement button.

SETTING	SPECIES		RESULT DISPLAY (readings are examples)
do	Dog	Canine	d 15
	Cat	Feline	
ho	Horse	Equine	H 18

REPLACING THE PROBE BASE

During the use some dirt may collect in the probe base, affecting the probe movement. The probe base should be replaced every 12 months or if the probe no longer moves smoothly or the probe base error (E 01 or E 03) is displayed.

1. Unscrew the probe base collar.
2. Remove the probe base by tilting the tonometer downward and pull the probe base out of the tonometer.
3. Put a new probe base into the tonometer.
4. Screw the collar in, to lock the probe base.

NOTE The probe base should be replaced with a new probe base every 12 months, or if the cleaning is not effective.



CLEANING THE PROBE BASE

The probe base should be cleaned every 6 months. The probe base can be reused after careful cleaning. Make sure the probe base is completely dry before using it.

1. Unscrew the probe base collar and remove the probe base from the tonometer.
2. Soak the probe base in 100% pure alcohol for 5 to 30 minutes.
3. Remove the probe base from alcohol and dry carefully by blowing clean canned or compressed air into the probe base. This will additionally remove possible residual dirt.
4. Insert the completely dry probe base and screw the collar back to the tonometer.

NOTE Do not use water or soap solution to clean the probe base.

CLEANING THE TONOMETER SURFACE

The tonometer may be surface-cleaned using a soft cloth dampened with either a commercial, nonabrasive cleaner or a solution of 70% alcohol in water.

Lightly wipe the surfaces of the tonometer and make sure surfaces are dry after cleaning, especially areas near the main buttons, the display and the probe cover.

Chemical resistance to the following liquids has been tested:

- Ethanol
- 2-propanol
- Mild soap solution
- 95% Pursept solution



CAUTION

Do not spray, pour or spill liquid onto the Icare tonometer, its accessories, connectors, switches or openings in the chassis.

Do not leave the surface of the tonometer wet; dry it with a soft cloth.
The probe holder/base has to be removed from the tonometer to be cleaned.