

# TONOVET *Pet*

icare  
TONOVET

## TONOMETER

The **TONOVET Pet** tonometer offers a cost-effective and user-friendly solution for monitoring intraocular pressure (IOP) in pets from the comfort of their home. Designed with reliable rebound technology, this device ensures effortless and accurate measurements without the need for anesthesia or calibration.




- Cost-effective option for IOP home monitoring by pet owner
- Reliable rebound technology for easy at home-use
- Results are displayed in 4 ranges

<15  
mmHg

15-25  
mmHg

26-35  
mmHg

>35  
mmHg

- The 4 measurement ranges allow the pet owner to cooperate with their veterinarian on their pet's treatment plan, providing a simplified way to monitor their pet's IOP without worrying about normal daily fluctuations.
- TONOVET Pet is designed for sale or lease to the pet owner by their veterinarian
- No topical anesthesia and no calibration
- Canine Setting only 

With the **TONOVET Pet**, pet owners can actively participate in managing their pet's eye health, working closely with veterinarians to ensure the best care possible.



NO ANESTHESIA  
NO CALIBRATION

[www.tonovet.com](http://www.tonovet.com)

